

兔抗 HMGN2 多克隆抗体

中文名称:	兔抗 HMGN2 多克隆抗体
英文名称:	Anti-HMGN2 rabbit polyclonal antibody
别名:	high mobility group nucleosomal binding domain 2; HMG17
相关类别	一抗
储存	冷冻 (-20℃)
宿主	Rabbit
抗原	HMGN2
反应种属	Human, Mouse, Rat
标记物	Unconjugate
克隆类型	rabbit polyclonal

技术规格

Background:	The protein encoded by this gene binds nucleosomal DNA and is associated with transcriptionally active chromatin. Along with a similar protein, HMGN1, the encoded protein may help maintain an open chromatin configuration around transcribable genes. The protein has also been found to have antimicrobial activity against bacteria, viruses and fungi.
Applications:	ELISA, IHC
Name of antibody:	HMGN2
Immunogen:	Fusion protein of human HMGN2
Full name:	high mobility group nucleosomal binding domain 2
Synonyms:	HMG17

SwissProt:	P05204
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommended dilution:	30-150

