

兔抗 HMGN5 多克隆抗体

中文名称：兔抗 HMGN5 多克隆抗体

英文名称：Anti-HMGN5 rabbit polyclonal antibody

别名：NSBP1, NBP-45

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HMGN5

反应种属：Human


标记物：Unconjugate



克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a nuclear protein with similarities to the high mobility group proteins, HMG14 and HMG17, which suggests that this protein may function as a nucleosomal binding and transcriptional activating protein.
Applications:	ELISA, IHC
Name of antibody:	HMGN5
Immunogen:	Synthetic peptide of human HMGN5
Full name:	High mobility group nucleosome binding domain 5
Synonyms :	NSBP1, NBP-45
SwissProt:	P82970

ELISA Recommended dilution:	1000-5000
IHC positive control:	Human brain
IHC Recommend dilution:	25-100





mlbio 酶联生物
Good elisakit producers