

兔抗 MAP3K6 多克隆抗体

中文名称：兔抗 MAP3K6 多克隆抗体

英文名称： Anti-MAP3K6 rabbit polyclonal antibody

别名： ASK2; MEKK6; MAPKKK6

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： MAP3K

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a serine/threonine protein kinase that forms a component of protein kinase-mediated signal transduction cascades. The encoded kinase participates in the regulation of vascular endothelial growth factor (VEGF) expression. Alternative splicing results in multiple transcript variants.
Applications:	IHC
Name of antibody:	MAP3K6
Immunogen:	Synthesized peptide derived from internal of human MAP3K6
Full name:	mitogen-activated protein kinase kinase kinase 6
Synonyms :	ASK2; MEKK6; MAPKKK6

SwissProt:	O95382
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100

