

兔抗 MEF2B 多克隆抗体

中文名称: 兔抗 MEF2B 多克隆抗体

英文名称: Anti-MEF2B rabbit polyclonal antibody

别 名: RSRFR2

相关类别: 一抗

储 存: 冷冻(-20℃) 避光

抗 原: MEF2B

宿 主: Rabbit

反应种属: Human Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

The product of this gene is a member of the MADS/MEF2 family of DNA binding proteins. The protein is thought to regulate gene expression, including expression of the smo oth muscle myosin heavy chain gene. This region undergo es considerable alternative splicing, with transcripts suppor ting two non-overlapping loci (GeneID 729991 and 100271 849) as well as numerous read-through transcripts that sp an both loci (annotated as GeneID 4207). Several isoforms of this protein are expressed from either this locus or from some of the read-through transcripts annotated on Gen eID 4207.



Applications:	WB
Name of antibody:	MEF2B
Immunogen:	Synthesized peptide derived from internal of human MEF2 B.
Full name:	myocyte enhancer factor 2B
Synonyms:	RSRFR2
SwissProt:	Q02080
WB Predicted band size:	39 kDa
WB Positive control:	HT-29 cells lysate
WB Recommended dilution:	500-3000

