

## 兔抗 MEOX2 多克隆抗体

中文名称：兔抗 MEOX2 多克隆抗体

英文名称： Anti-MEOX2 rabbit polyclonal antibody

别名： GAX; MOX2

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： MEOX2

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a member of a subfamily of non-clustered, diverged, antennapedia-like homeobox-containing genes. The encoded protein may play a role in the regulation of vertebrate limb myogenesis. Mutations in the related mouse protein may be associated with craniofacial and/or skeletal abnormalities, in addition to neurovascular dysfunction observed in Alzheimer's disease.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	MEOX2
<b>Immunogen:</b>	Synthesized peptide derived from internal of human MEOX2.
<b>Full name:</b>	mesenchyme homeobox 2

<b>Synonyms :</b>	GAX; MOX2
<b>SwissProt:</b>	P50222
<b>WB Predicted band size:</b>	34 kDa
<b>WB Positive control:</b>	COLO205 cells lysate
<b>WB Recommended dilution</b> :	500-3000

