

兔抗 MERTK 多克隆抗体

中文名称：兔抗 MERTK 多克隆抗体

英文名称： Anti-MERTK rabbit polyclonal antibody

抗 原： MERTK

别 名： MER, RP38, c-mer

储 存： 冷冻（-20℃）

宿 主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a member of the MER/AXL/TYRO3 receptor kinase family and encodes a transmembrane protein with two fibronectin type-III domains, two Ig-like C2-type (immunoglobulin-like) domains, and one tyrosine kinase domain. Mutations in this gene have been associated with disruption of the retinal pigment epithelium (RPE) phagocytosis pathway and onset of autosomal recessive retinitis pigmentosa (RP).
Applications:	ELISA, IHC
Name of antibody:	MERTK
Immunogen:	Fusion protein of human MERTK

Full name:	c-mer proto-oncogene tyrosine kinase
Synonyms :	MER, RP38, c-mer
SwissProt:	Q12866
ELISA Recommended dilution :	1000-5000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	25-100

