

兔抗 MEST 多克隆抗体

中文名称：兔抗 MEST 多克隆抗体

英文名称：Anti-MEST rabbit polyclonal antibody

别名：mesoderm specific transcript; PEG1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MEST

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of the alpha/beta hydrolase superfamily. It is imprinted, exhibiting preferential expression from the paternal allele in fetal tissues, and isoform-specific imprinting in lymphocytes. The loss of imprinting of this gene has been linked to certain types of cancer and may be due to promotor switching. The encoded protein may play a role in development. Alternatively spliced transcript variants encoding multiple isoforms have been identified for this gene. Pseudogenes of this gene are located on the short arm of chromosomes 3 and 4, and the long arm of chromosomes 6 and 15.

Applications:	ELISA, IHC
Name of antibody:	MEST
Immunogen:	Fusion protein of human MEST
Full name:	mesoderm specific transcript
Synonyms:	PEG1
SwissProt:	Q5EB52
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300

