

兔抗 MET(Phospho-Tyr1003) 多克隆抗体

中文名称：兔抗 MET(Phospho-Tyr1003) 多克隆抗体

英文名称： Anti-MET(Phospho-Tyr1003) rabbit polyclonal antibody

别名： HGFR; AUTS9; RCCP2; c-Met

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： MET(Phospho-Tyr1003)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The proto-oncogene MET product is the hepatocyte growth factor receptor and encodes tyrosine-kinase activity. The primary single chain precursor protein is post-translationally cleaved to produce the alpha and beta subunits, which are disulfide linked to form the mature receptor. Various mutations in the MET gene are associated with papillary renal carcinoma. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	WB
Name of antibody:	MET(Phospho-Tyr1003)
Immunogen:	Peptide sequence around phosphorylation site of tyrosine

	1003 (V-D-Y(p)-R-A) derived from Human c-Met.
Full name:	MET proto-oncogene, receptor tyrosine kinase
Synonyms :	HGFR; AUTS9; RCCP2; c-Met
SwissProt:	P08581
WB Predicted band size:	156 kDa
WB Positive control:	HepG2 cells lysates
WB Recommended dilution:	500-1000

