

兔抗 METTL17 多克隆抗体

中文名称: 兔抗 METTL17 多克隆抗体

英文名称: Anti-METTL17 rabbit polyclonal antibody

别名: methyltransferase like 17; METT11D1

抗原: METTL17

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	May be a component of the mitochondrial small ribosomal subunit.
Applications:	ELISA, WB, IHC
Name of antibody:	METTL17
Immunogen:	Synthetic peptide of human METTL17
Full name:	methyltransferase like 17
Synonyms:	METT11D1
SwissProt:	Q9H7H0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300

WB Predicted band size:	51 kDa
WB Positive control:	Human colorectal cancer tissue lysate
WB Recommended dilution:	500-2000



