

## 兔抗 METTL1 多克隆抗体

中文名称：兔抗 METTL1 多克隆抗体

英文名称：Anti-METTL1 rabbit polyclonal antibody

别名：methyltransferase like 1; TRM8; TRMT8; C12orf1; YDL201w

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：METTL1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is similar in sequence to the <i>S. cerevisiae</i> YDL201w gene. The gene product contains a conserved S-adenosylmethionine-binding motif and is inactivated by phosphorylation. Alternative splice variants encoding different protein isoforms have been described for this gene. A pseudogene has been identified on chromosome X.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	METTL1
<b>Immunogen:</b>	Fusion protein of human METTL1
<b>Full name:</b>	methyltransferase like 1
<b>Synonyms:</b>	TRM8; TRMT8; C12orf1; YDL201w

<b>SwissProt:</b>	Q9UBP6
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	50-300



