

兔抗 MFAP3L 多克隆抗体

中文名称：兔抗 MFAP3L 多克隆抗体

英文名称： Anti-MFAP3L rabbit polyclonal antibody

别名： NYD-sp9

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： MFAP3L

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

MFAP3L (microfibrillar-associated protein 3-like), also known as HSD39 or testis development protein NYD-SP9, is a 409 amino acid single-pass type I cell membrane protein that contains one Ig-like (immunoglobulin-like) domain. Found primarily in testis, MFAP3L is encoded by a gene that is located on chromosome 4 and is expressed as three isoforms due to alternative splicing events. Representing approximately 6% of the human genome, chromosome 4 contains nearly 900 genes, one of which is the Huntington gene, which is found to encode an expanded glutamine tract in cases of Huntington's disease.

Applications:	ELISA, IHC
Name of antibody:	MFAP3L
Immunogen:	Fusion protein of human MFAP3L
Full name:	microfibrillar-associated protein 3-like
Synonyms :	NYD-sp9
SwissProt:	O75121
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human colon cancer and human liver cancer
IHC Recommend dilution:	100-300

