

兔抗 MGST1 多克隆抗体

中文名称：兔抗 MGST1 多克隆抗体

英文名称：Anti-MGST1 rabbit polyclonal antibody

别名：MGST; GST12; MGST-I

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MGST1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The MAPEG (Membrane Associated Proteins in Eicosanoid and Glutathione metabolism) family consists of six human proteins, two of which are involved in the production of leukotrienes and prostaglandin E, important mediators of inflammation. Other family members, demonstrating glutathione S-transferase and peroxidase activities, are involved in cellular defense against toxic, carcinogenic, and pharmacologically active electrophilic compounds.
Applications:	ELISA, IHC
Name of antibody:	MGST1

Immunogen:	Synthetic peptide of human MGST1
Full name:	microsomal glutathione S-transferase 1
Synonyms :	MGST; GST12; MGST-I
SwissProt:	P10620
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	25-100

