

## 兔抗 MGST2 多克隆抗体

中文名称：兔抗 MGST2 多克隆抗体

英文名称： Anti-MGST2 rabbit polyclonal antibody

别名： GST2; MGST-II

抗原： MGST2

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

|                          |  |
|--------------------------|--|
| <b>Background:</b>       | The MAPEG (Membrane Associated Proteins in Eicosanoid and Glutathione metabolism) family consists of six human proteins, several of which are involved in the production of leukotrienes and prostaglandin E, important mediators of inflammation. This gene encodes a protein which catalyzes the conjugation of leukotriene A4 and reduced glutathione to produce leukotriene C4. Alternatively spliced transcript variants encoding different isoforms have been identified in this gene. |
| <b>Applications:</b>     | WB   |
| <b>Name of antibody:</b> | MGST2  |

|                                 |  |
|---------------------------------|--|
| <b>Immunogen:</b>               | Synthesized peptide derived from Internal of human MGST 2. |
| <b>Full name:</b>               | microsomal glutathione S-transferase 2                     |
| <b>Synonyms :</b>               | GST2; MGST-II  |
| <b>SwissProt:</b>               | Q99735   |
| <b>WB Predicted band size:</b>  | 17 kDa   |
| <b>WB Positive control:</b>     | COLO205 cells lysate                                       |
| <b>WB Recommended dilution:</b> | 500-3000   |

