

## 兔抗 MIA 多克隆抗体

中文名称: 兔抗 MIA 多克隆抗体

英文名称: Anti-MIA rabbit polyclonal antibody

别名: melanoma inhibitory activity; CD-RAP

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: MIA

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MIA
<b>Immunogen:</b>	Fusion protein of human MIA
<b>Full name:</b>	melanoma inhibitory activity
<b>Synonyms:</b>	CD-RAP
<b>SwissProt:</b>	Q16674
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer and Human colorectal cancer

IHC Recommend dilution: 100-300



