

兔抗 MICU3 多克隆抗体

中文名称：兔抗 MICU3 多克隆抗体

英文名称：Anti-MICU3 rabbit polyclonal antibody

别名：EFHA2

相关类别：一抗

抗原：MICU3

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May play a role in mitochondrial calcium uptake.
Applications:	ELISA, IHC
Name of antibody:	MICU3
Immunogen:	Synthetic peptide of human MICU3
Full name:	mitochondrial calcium uptake family member 3
Synonyms:	EFHA2
SwissProt:	Q86XE3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer;Human prostate cancer
IHC Recommend dilution:	25-100

