

## 兔抗 MIEN1 多克隆抗体

中文名称：兔抗 MIEN1 多克隆抗体

英文名称：Anti-MIEN1 rabbit polyclonal antibody

级别：C35; ORB3; XTP4; RDX12; C17orf37

抗原：MIEN1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

相关类别：一抗

### 技术规格

<b>Background:</b>	Increases cell migration by inducing filopodia formation at the leading edge of migrating cells. Plays a role in regulation of apoptosis, possibly through control of CASP3. May be involved in a redox-related process.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MIEN1
<b>Immunogen:</b>	Synthetic peptide of human MIEN1
<b>Full name:</b>	migration and invasion enhancer 1
<b>Synonyms:</b>	C35; ORB3; XTP4; RDX12; C17orf37
<b>SwissProt:</b>	Q9BRT3
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human prostate cancer
<b>IHC Recommend dilution:</b>	20-100



