

兔抗 MIER2 多克隆抗体

中文名称：兔抗 MIER2 多克隆抗体

英文名称： Anti-MIER2 rabbit polyclonal antibody

别名： Mi-er2; KIAA1193

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MIER2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	MIER2(mesoderm induction early response 1, family member 2), also known as Mi-er2 or KIAA1193, is a protein that in humans is encoded by the MIER2 gene. It is a transcriptional repressor.
Applications:	ELISA, WB, IHC
Name of antibody:	MIER2
Immunogen:	Fusion protein of human MIER2
Full name:	mesoderm induction early response 1, family member 2
Synonyms :	Mi-er2; KIAA1193
SwissProt:	Q8N344
ELISA Recommended dilution:	1000-2000

IHC positive control:	Human colon cancer and human liver cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	60 kDa
WB Positive control:	Hela cells
WB Recommended dilution:	200-1000



