

兔抗 MIIP 多克隆抗体

中文名称: 兔抗 MIIP 多克隆抗体

英文名称: Anti-MIIP rabbit polyclonal antibody

别名: IIP45

储存: 冷冻 (-20℃)

抗原: MIIP

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a protein that interacts with the onco gene protein insulin-like growth factor binding protein 2 and may function as an inhibitor of cell migration and invasion. This protein also interacts with the cell division protein 20 and may be involved in regulating mitotic progression. This protein may function as a tumor suppressor by inhibiting the growth of certain cancers.
Applications:	ELISA, IHC
Name of antibody:	MIIP
Immunogen:	Synthetic peptide of human MIIP

Full name:	migration and invasion inhibitory protein
Synonyms:	IP45
SwissProt:	Q5JXC2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human esophagus cancer
IHC Recommend dilution:	25-100



