

兔抗 MIP 多克隆抗体

中文名称：兔抗 MIP 多克隆抗体

英文名称：Anti-MIP rabbit polyclonal antibody

别名：AQP0; LIM1; MP26; MIP26; CTRCT15

相关类别：一抗

抗原：MIP

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Major intrinsic protein is a member of the water-transporting aquaporins as well as the original member of the MIP family of channel proteins. The function of the fiber cell membrane protein encoded by this gene is undetermined, yet this protein is speculated to play a role in intracellular communication. The MIP protein is expressed in the ocular lens and is required for correct lens function. This gene has been mapped among aquaporins AQP2, AQP5, and AQP6, in a potential gene cluster at 12q13.
Applications:	ELISA, IHC
Name of antibody:	MIP

Immunogen:	Synthetic peptide of humanMIP
Full name:	major intrinsic protein of lens fiber
Synonyms :	AQP0; LIM1; MP26; MIP26; CTRCT15
SwissProt:	P30301
ELISA Recommended dilution:	500-5000
IHC positive control:	Human tonsil and Human thyroid cancer
IHC Recommend dilution:	50-200

