

兔抗 MKKS 多克隆抗体

中文名称：兔抗 MKKS 多克隆抗体

英文名称： Anti-MKKS rabbit polyclonal antibody

别名： KMS; MKS; BBS6; HMCS

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MKKS

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a protein which shares sequence similarity with other members of the type II chaperonin family . The encoded protein is a centrosome-shuttling protein and plays an important role in cytokinesis. This protein also interacts with other type II chaperonin members to form a complex known as the BBSome, which involves ciliary membrane biogenesis. This protein is encoded by a downstream open reading frame (dORF). Several upstream open reading frames (uORFs) have been identified, which repress the translation of the dORF, and two of which can encode small mitochondrial membrane proteins.

Applications:	ELISA, IHC
Name of antibody:	MKKS
Immunogen:	Fusion protein of human MKKS
Full name:	McKusick-Kaufman syndrome
Synonyms :	KMS; MKS; BBS6; HMCS
SwissProt:	Q9NPJ1
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer and human brain
IHC Recommend dilution:	50-200

