

## 兔抗 MKNK1 多克隆抗体

中文名称：兔抗 MKNK1 多克隆抗体

英文名称： Anti-MKNK1 rabbit polyclonal antibody

别名： MNK1

抗原： MKNK1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	May play a role in the response to environmental stress and cytokines. Appears to regulate translation by phosphorylating EIF4E, thus increasing the affinity of this protein for the 7-methylguanosine-containing mRNA cap.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	MKNK1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human MKNK1.
<b>Full name:</b>	MAP kinase interacting serine/threonine kinase 1
<b>Synonyms :</b>	MNK1
<b>SwissProt:</b>	Q9BUB5

<b>WB Predicted band size:</b>	51 kDa
<b>WB Positive control:</b>	Huvec cells lysate
<b>WB Recommended dilution:</b>	500-3000

