

兔抗 MLF1 多克隆抗体

中文名称：兔抗 MLF1 多克隆抗体

英文名称： Anti-MLF1 rabbit polyclonal antibody

别名： myeloid leukemia factor 1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MLF1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an oncoprotein which is thought to play a role in the phenotypic determination of hemopoietic cells. Translocations between this gene and nucleophosmin have been associated with myelodysplastic syndrome and acute myeloid leukemia. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB
Name of antibody:	MLF1
Immunogen:	Fusion protein of human MLF1
Full name:	myeloid leukemia factor 1

SwissProt:	P58340
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	31 kDa
WB Positive control:	Mouse heart tissue, Human cerebrum tissue, Human cerebella tissue lysates
WB Recommended dilution:	500-2000

