

兔抗 MLNR 多克隆抗体

中文名称：兔抗 MLNR 多克隆抗体

英文名称：Anti-MLNR rabbit polyclonal antibody

别名：GPR38; MTLR1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MLNR

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Motilin is a 22 amino acid peptide hormone expressed throughout the gastrointestinal (GI) tract. The protein encoded by this gene is a motilin receptor which is a member of the G-protein coupled receptor 1 family. This member is a multi-pass transmembrane protein, and is an important therapeutic target for the treatment of hypomotility disorders.
Applications:	ELISA, IHC
Name of antibody:	MLNR
Immunogen:	Synthetic peptide of human MLNR
Full name:	motilin receptor

Synonyms :	GPR38; MTLR1
SwissProt:	O43193
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

