

兔抗 MLN 多克隆抗体

中文名称:	兔抗 MLN 多克隆抗体
英文名称:	Anti-MLN rabbit polyclonal antibody
别名:	motilin
相关类别:	一抗
储存:	冷冻 (-20℃)
宿主:	Rabbit
抗原:	MLN
反应种属:	Human
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	This gene encodes a small peptide hormone that is secreted by cells of the small intestine to regulate gastrointestinal contractions and motility. Proteolytic processing of the secreted protein produces the mature peptide and a byproduct referred to as motilin-associated peptide (MAP). Three transcript variants encoding different preproprotein isoforms but the same mature peptide have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	MLN
Immunogen:	Fusion protein of human MLN

Full name:	motilin
SwissProt:	P12872
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-100

