

EFNB2 抗原（重组蛋白）

中文名称：EFNB2 抗原（重组蛋白）

英文名称：EFNB2 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 28-229 amino acids of human EFNB2

技术规格

Full name:	ephrin-B2
Synonyms:	HTKL, EPLG5, Htk-L, LERK5
Swissprot:	P52799
Gene Accession:	BC069342
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the ephrin (EPH) family. The ephrins and EPH-related receptors comprise the largest subfamily of receptor protein-tyrosine kinases and have been implicated in mediating developmental events, especially in the nervous system and in erythropoiesis. Based on their structures and sequence relationships, ephrins are divided into the ephrin-A (EFNA) class, which are anchored to the membrane by a glycosylphosphatidylinositol linkage, and the ephrin-B (EFNB) class, which are transmembrane proteins. This gene encodes an EFNB class ephrin which binds to the EPHB4 and EPHA3 receptors.

mlbio 梅联生物
Good elisakit producers