

EGFL7 抗原（重组蛋白）

中文名称：EGFL7 抗原（重组蛋白）

英文名称：EGFL7 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 24-273 amino acids of human EGFL7

技术规格

Full name:	EGF-like-domain, multiple 7
Synonyms:	NEU1, ZNEU1, VE-STATIN, RP11-251M1.2
Swissprot:	Q9UHF1
Gene Accession:	BC012377
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a secreted endothelial cell protein that contains two epidermal growth factor-like domains. The encoded protein may play a role in regulating vasculogenesis. This protein may be involved in the growth and proliferation of tumor cells. Alternate splicing results in multiple transcript variants.