

EGR3 抗原（重组蛋白）

中文名称：EGR3 抗原（重组蛋白）

英文名称：EGR3 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to C terminal 300 amino acids of human EGR3

技术规格

Full name:	early growth response 3
Synonyms:	EGR-3; PILOT
Swissprot:	Q06889
Gene Accession:	BC112279
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a transcriptional regulator that belongs to the EGR family of C2H2-type zinc-finger proteins. It is an immediate-early growth response gene which is induced by mitogenic stimulation. The protein encoded by this gene participates in the transcriptional regulation of genes in controlling biological rhythm. It may also play a role in a wide variety of processes including muscle development, lymphocyte development, endothelial cell growth and migration, and neuronal development.