

ETS1 抗原（重组蛋白）

中文名称：ETS1 抗原（重组蛋白）

英文名称：ETS1 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to C terminal 209 amino acids of human ETS1

技术规格

Full name:	v-ets erythroblastosis virus E26 oncogene homolog 1 (avian)
Synonyms:	ETS-1, EWSR2
Swissprot:	P14921
Gene Accession:	BC017314
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the ETS family of transcription factors, which are defined by the presence of a conserved ETS DNA-binding domain that recognizes the core consensus DNA sequence GGAA/T in target genes. These proteins function either as transcriptional activators or repressors of numerous genes, and are involved in stem cell development, cell senescence and death, and tumorigenesis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.