

FABP4 抗原（重组蛋白）

中文名称：FABP4 抗原（重组蛋白）

英文名称：FABP4 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 1-132 amino acids of human FABP4

技术规格

Full name:	fatty acid binding protein 4
Synonyms:	aP2; ALBP; AFABP; A-FABP; HEL-S-104
Swissprot:	P15090
Gene Accession:	BC003672
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	FABP4 encodes the fatty acid binding protein found in adipocytes. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. It is thought that FABPs roles include fatty acid uptake, transport, and metabolism.