

## FOS 抗原（重组蛋白）

中文名称：FOS 抗原（重组蛋白）

英文名称： FOS Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 131-380 amino acids of human FOS

技术规格

<b>Full name:</b>	Fos proto-oncogene, AP-1 transcription factor subunit
<b>Synonyms:</b>	p55; AP-1; C-FOS
<b>Swissprot:</b>	P01100
<b>Gene Accession:</b>	BC004490
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differen