

FOXS1 抗原（重组蛋白）

中文名称：FOXS1 抗原（重组蛋白）

英文名称：FOXS1 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 131-330 amino acids of human FOXS1

技术规格

Full name:	forkhead box S1
Synonyms:	FKHL18, FREAC10
Swissprot:	O43638
Gene Accession:	BC013408
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The forkhead family of transcription factors belongs to the winged helix class of DNA-binding proteins. The protein encoded by this intronless gene contains a forkhead domain and is found predominantly in aorta and kidney. The function of the encoded protein is unknown.