

## GAS7 抗原（重组蛋白）

中文名称：GAS7 抗原（重组蛋白）

英文名称：GAS7 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 184-382 amino acids of human GAS7

技术规格

<b>Full name:</b>	growth arrest-specific 7
<b>Synonyms:</b>	MLL/GAS7
<b>Swissprot:</b>	O60861
<b>Gene Accession:</b>	BC001152
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Growth arrest-specific 7 is expressed primarily in terminally differentiated brain cells and predominantly in mature cerebellar Purkinje neurons. GAS7 plays a putative role in neuronal development. Several transcript variants encoding proteins which vary in the N-terminus have been described.