

兔抗 GSR 多克隆抗体

中文名称：兔抗 GSR 多克隆抗体

英文名称： Anti-GSR rabbit polyclonal antibody

抗 原： GSR

储 存： 冷冻（-20℃）

宿 主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the class-I pyridine nucleotide-disulfide oxidoreductase family. This enzyme is a homodimeric flavoprotein. It is a central enzyme of cellular antioxidant defense, and reduces oxidized glutathione disulfide (GSSG) to the sulfhydryl form GSH, which is an important cellular antioxidant. Rare mutations in this gene result in hereditary glutathione reductase deficiency. Multiple alternatively spliced transcript variants encoding different isoforms have been found.
Applications:	ELISA, WB, IHC
Name of antibody:	GSR
Immunogen:	Fusion protein of human GSR
Full name:	Glutathione reductase

SwissProt:	P00390
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human cervical cancer and Human ovarian cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	56 kDa
WB Positive control:	Jurkat, Hela and A549 cells
WB Recommended dilution:	1000-5000



