

兔抗 HNRNPCL1 多克隆抗体

- 中文名称：兔抗 HNRNPCL1 多克隆抗体
- 英文名称：Anti-HNRNPCL1 rabbit polyclonal antibody
- 别名：heterogeneous nuclear ribonucleoprotein C-like 1; HNRNPCL1
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：HNRNPCL1
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	May play a role in nucleosome assembly by neutralizing basic proteins such as A and B core hnRNPs.
Applications:	ELISA, IHC
Name of antibody:	HNRNPCL1
Immunogen:	Synthetic peptide of human HNRNPCL1
Full name:	heterogeneous nuclear ribonucleoprotein C-like 1
Synonyms:	HNRNPCL1
SwissProt:	O60812
ELISA Recommended dilution:	5000-10000
IHC positive control	Human esophagus cancer and Human tonsil

:	
IHC Recommend dilution:	25-100



