

兔抗 GSS 多克隆抗体

中文名称：兔抗 GSS 多克隆抗体

英文名称：Anti-GSS rabbit polyclonal antibody

别名：GSS; GSXS; MGC14098

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：GSS

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Glutathione is important for a variety of biological functions, including protection of cells from oxidative damage by free radicals, detoxification of xenobiotics, and membrane transport. The protein encoded by this gene functions as a homodimer to catalyze the second step of glutathione biosynthesis, which is the ATP-dependent conversion of gamma-L-glutamyl-L-cysteine to glutathione. Defects in this gene are a cause of glutathione synthetase deficiency.
Applications:	WB
Name of antibody:	GSS

Immunogen:	Fusion protein of human GSS
Full name:	glutathione synthetase
Synonyms :	GSS; GSHS; MGC14098
SwissProt:	P48637
WB Predicted band size:	52 kDa
WB Positive control:	293T, HeLa and HT-1080 cells
WB Recommended dilution:	500-2000

