

## GNAZ 抗原（重组蛋白）

中文名称：GNAZ 抗原（重组蛋白）

英文名称：GNAZ Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 106-355 amino acids of human GNAZ

技术规格

<b>Full name:</b>	G protein subunit alpha z
<b>Synonyms:</b>	gz-alpha
<b>Swissprot:</b>	P19086
<b>Gene Accession:</b>	BC026342
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is a member of a G protein subfamily that mediates signal transduction in pertussis toxin-insensitive systems. This encoded protein may play a role in maintaining the ionic balance of perilymphatic and endolymphatic cochlear fluids.