

兔抗 GSTA2 多克隆抗体

中文名称：兔抗 GSTA2 多克隆抗体

英文名称： Anti-GSTA2 rabbit polyclonal antibody

别名： GST2; GTA2; GTH2; GSTA2-2

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： GSTA2

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Cytosolic and membrane-bound forms of glutathione S-transferase are encoded by two distinct supergene families. These enzymes function in the detoxification of electrophilic compounds, including carcinogens, therapeutic drugs, environmental toxins and products of oxidative stress, by conjugation with glutathione. The genes encoding these enzymes are known to be highly polymorphic. These genetic variations can change an individual's susceptibility to carcinogens and toxins as well as affect the toxicity and efficacy of some drugs.
Applications:	ELISA, IHC

Name of antibody:	GSTA2
Immunogen:	Fusion protein of human GSTA2
Full name:	glutathione S-transferase alpha 2
Synonyms :	GST2; GTA2; GTH2; GSTA2-2
SwissProt:	P09210
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human brain and human lung cancer
IHC Recommend dilution:	50-200

