

GPC6 抗原（重组蛋白）

中文名称：GPC6 抗原（重组蛋白）

英文名称：GPC6 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 229-529 amino acids of human GPC6

技术规格

Full name:	glypican 6
Synonyms:	OMIMD1
Swissprot:	Q9Y625
Gene Accession:	BC106947
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The glypicans comprise a family of glycosylphosphatidylinositol-anchored heparan sulfate proteoglycans, and they have been implicated in the control of cell growth and cell division. The glypican encoded by this gene is a putative cell surface coreceptor for growth factors, extracellular matrix proteins, proteases and anti-proteases.