

兔抗 HNRNPD(Phospho-Ser83) 多克隆抗体

中文名称：兔抗 HNRNPD(Phospho-Ser83) 多克隆抗体

英文名称：Anti-HNRNPD(Phospho-Ser83) rabbit polyclonal antibody

别名：P37; AUF1; AUF1A; HNRPD; hnRNPD0

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：HNRNPD(Phospho-Ser83)

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are nuclear acid binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm.
Applications:	WB, IHC, IF

Name of antibody:	HNRNPD(Phospho-Ser83)
Immunogen:	Peptide sequence around phosphorylation site of Serine 83(N-S-S(p)-P-R) derived from Human hnRPD .
Full name:	heterogeneous nuclear ribonucleoprotein D (AU-rich element RNA binding protein 1, 37kDa)
Synonyms :	P37; AUF1; AUF1A; HNRPD; hnRNP0
SwissProt:	Q14103
IHC positive control :	Human lung carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	38 kDa
WB Positive control :	Jurkat cells lysates
WB Recommended dilution:	500-1000
IF positive control:	Hela cells
IF Recommend dilution:	100-200



