

兔抗 MNX1 多克隆抗体

中文名称：兔抗 MNX1 多克隆抗体

英文名称：Anti-MNX1 rabbit polyclonal antibody

别名：motor neuron and pancreas homeobox 1; HB9; HLXB9; SCRA1; HOXHB9

相关类别：一抗

抗原：MNX1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a nuclear protein, which contains a homeobox domain and is a transcription factor. Mutations in this gene result in Currarino syndrome, an autosomal dominant congenital malformation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB
Name of antibody:	MNX1
Immunogen:	Synthetic peptide of human MNX1
Full name:	motor neuron and pancreas homeobox 1
Synonyms:	HB9; HLXB9; SCRA1; HOXHB9

SwissProt:	P50219
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	41 kDa
WB Positive control:	Raji and HepG2 cell lysates
WB Recommended dilution:	200-1000

