

## 兔抗 HNRNPH1 多克隆抗体


- 中文名称：兔抗 HNRNPH1 多克隆抗体
- 英文名称：Anti-HNRNPH1 rabbit polyclonal antibody
- 别名：HNRPH; HNRPH1; hnRNPH
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：HNRNPH1
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

**Background:**

This gene encodes a member of a subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins that complex with heterogeneous nuclear RNA. These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some may shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has three repeats of quasi-RRM domains that bind to RNA and is very similar to the family member HNRPF. This

	gene may be associated with hereditary lymphedema type I.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	HNRNPH1
<b>Immunogen:</b>	Synthetic peptide of human HNRNPH1
<b>Full name:</b>	heterogeneous nuclear ribonucleoprotein H1 (H)
<b>Synonyms :</b>	HNRPH; HNRPH1; hnRNPH
<b>SwissProt:</b>	P31943
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size :</b>	49 kDa
<b>WB Positive control:</b>	LoVo and A375 cells, mouse brain tissue, HepG2, A431 and HeLa cells
<b>WB Recommended dilution:</b>	200-1000



**mlbio** 梅联生物  
Good elisakit producers