

兔抗 GADD45A 多克隆抗体

中文名称：兔抗 GADD45A 多克隆抗体

英文名称： Anti-GADD45A rabbit polyclonal antibody

别名： DDIT1; GADD45

相关类别： 一抗

抗原： GADD45A

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. The protein encoded by this gene responds to environmental stresses by mediating activation of the p38/JNK pathway via MTK1/MEKK4 kinase. The DNA damage-induced transcription of this gene is mediated by both p53-dependent and -independent mechanisms. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
Applications:	ELISA, IHC

Name of antibody:	GADD45A
Immunogen:	Fusion protein of human GADD45A
Full name:	growth arrest and DNA-damage-inducible, alpha
Synonyms :	DDIT1; GADD45
SwissProt:	P24522
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human brain and human cervical cancer
IHC Recommend dilution:	25-100

