

兔抗 MOB3B 多克隆抗体

中文名称: 兔抗 MOB3B 多克隆抗体

英文名称: Anti-MOB3B rabbit polyclonal antibody

别名: MOB kinase activator 3B; MOB1D; C9orf35; MOBKL2B

抗原: MOB3B

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene shares similarity with the yeast Mob1 protein. Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation. This gene is located on the opposite strand as the interferon kappa precursor (IFNK) gene.
Applications:	ELISA, IHC
Name of antibody:	MOB3B
Immunogen:	Synthetic peptide of human MOB3B
Full name:	MOB kinase activator 3B
Synonyms:	MOB1D; C9orf35; MOBKL2B

SwissProt:	Q86TA1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	40-200

