

兔抗 GSTO1 多克隆抗体

中文名称：兔抗 GSTO1 多克隆抗体

英文名称： Anti-GSTO1 rabbit polyclonal antibody

别名： P28; SPG-R; GSTO 1-1; GSTTLp28; HEL-S-21

相关类别： 一抗

抗原： GSTO1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is an omega class glutathione S-transferase (GST) with glutathione-dependent thiol transferase and dehydroascorbate reductase activities. GSTs are involved in the metabolism of xenobiotics and carcinogens. The encoded protein acts as a homodimer and is found in the cytoplasm. Three transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	GSTO1
Immunogen:	Fusion protein of human GSTO1
Full name:	glutathione S-transferase omega 1

Synonyms :	P28; SPG-R; GSTO 1-1; GSTTLp28; HEL-S-21
SwissProt:	P78417
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human colon cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	28 kDa
WB Positive control:	Human liver cancer tissue
WB Recommended dilution:	500-2000



