

## 兔抗 MOB3C 多克隆抗体

中文名称：兔抗 MOB3C 多克隆抗体

英文名称： Anti-MOB3C rabbit polyclonal antibody

别名： MOB1E; MOBKL2C

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： MOB3C

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is similar to the yeast Mob1 protein. Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation. Alternatively spliced transcript variants encoding distinct isoforms have been observed.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	MOB3C
<b>Immunogen:</b>	Synthesized peptide derived from internal of human MOB3C.
<b>Full name:</b>	MOB kinase activator 3C
<b>Synonyms :</b>	MOB1E; MOBKL2C
<b>SwissProt:</b>	Q70IA8

<b>WB Predicted band size:</b>	26 kDa
<b>WB Positive control:</b>	RAW264.7 cells and K562 cells lysates
<b>WB Recommended dilution:</b>	500-3000

