

## HDAC10 抗原（重组蛋白）

中文名称：HDAC10 抗原（重组蛋白）

英文名称：HDAC10 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to N terminal 300 amino acids of human HDAC10

技术规格

<b>Full name:</b>	histone deacetylase 10
<b>Synonyms:</b>	HD10
<b>Swissprot:</b>	Q969S8
<b>Gene Accession:</b>	BC125083
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene belongs to the histone deacetylase family, members of which deacetylate lysine residues on the N-terminal part of the core histones. Histone deacetylation modulates chromatin structure, and plays an important role in transcriptional regulation, cell cycle progression, and developmental events. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.