

兔抗 MOB4 多克隆抗体

中文名称：兔抗 MOB4 多克隆抗体

英文名称：Anti-MOB4 rabbit polyclonal antibody

别名：2C4D; MOB1; MOB3; PHOCN; PREI3; CGI-95; MOBKL3

相关类别：一抗

抗原：MOB4

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene was identified based on its similarity with the mouse counterpart. Studies of the mouse counterpart suggest that the expression of this gene may be regulated during oocyte maturation and preimplantation following zygotic gene activation. Alternatively spliced transcript variants encoding distinct isoforms have been observed. Naturally occurring read-through transcription occurs between this locus and the neighboring locus HSPE1.
Applications:	ELISA, WB, IHC
Name of antibody:	MOB4
Immunogen:	Full length fusion protein

Full name:	MOB family member 4, phocein
Synonyms :	2C4D; MOB1; MOB3; PHOCN; PREI3; CGI-95; MOBKL3
SwissProt:	Q9Y3A3
ELISA Recommended dilution :	5000-10000
IHC positive control:	Human lung cancer and human breast cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	26 kDa
WB Positive control:	Human fetal brain tissue , Human cerebrum tissue , Mouse brain tissue
WB Recommended dilution:	500-2000



