

## 兔抗 MOCS3 多克隆抗体

中文名称：兔抗 MOCS3 多克隆抗体

英文名称： Anti-MOCS3 rabbit polyclonal antibody

别名： molybdenum cofactor synthesis 3; UBA4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MOCS3

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Molybdenum cofactor (MoCo) is necessary for the function of all molybdoenzymes. The protein encoded by this gene adenylates and activates molybdopterin synthase, an enzyme required for biosynthesis of MoCo. This gene contains no introns. A pseudogene of this gene is present on chromosome 14.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	MOCS3
<b>Immunogen:</b>	Fusion protein of human MOCS3
<b>Full name:</b>	molybdenum cofactor synthesis 3
<b>Synonyms:</b>	UBA4

<b>SwissProt:</b>	O95396
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	50 kDa
<b>WB Positive control:</b>	Mouse heart tissue lysate
<b>WB Recommended dilution:</b>	200-1000

