

兔抗 HNRNPUL1 多克隆抗体

- 中文名称: 兔抗 HNRNPUL1 多克隆抗体
- 英文名称: Anti-HNRNPUL1 rabbit polyclonal antibody
- 别名: heterogeneous nuclear ribonucleoprotein U like 1; E1BAP5; E1B-AP5; HNRNPUL1
- 相关类别: 一抗
- 储存: 冷冻 (-20℃)
- 宿主: Rabbit
- 抗原: HNRNPUL1
- 反应种属: Human, Mouse
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a nuclear RNA-binding protein of the heterogeneous nuclear ribonucleoprotein (hnRNP) family. This protein binds specifically to adenovirus early-1B-55kDa oncoprotein. It may play an important role in nucleocytoplasmic RNA transport, and its function is modulated by early-1B-55kDa in adenovirus-infected cells.
Applications:	ELISA, IHC
Name of antibody:	HNRNPUL1
Immunogen:	Synthetic peptide of human HNRNPUL1
Full name:	heterogeneous nuclear ribonucleoprotein U like 1

Synonyms:	E1BAP5; E1B-AP5; HNRPUL1
SwissProt:	Q9BUJ2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution :	20-100

