

## 兔抗 GAK 多克隆抗体

中文名称：兔抗 GAK 多克隆抗体

英文名称： Anti-GAK rabbit polyclonal antibody

别名： DNAJ26; DNAJC26

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： GAK

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	In all eukaryotes, the cell cycle is governed by cyclin-dependent protein kinases (CDKs), whose activities are regulated by cyclins and CDK inhibitors in a diverse array of mechanisms that involve the control of phosphorylation and dephosphorylation of Ser, Thr or Tyr residues. Cyclins are molecules that possess a consensus domain called the 'cyclin box.' In mammalian cells, 9 cyclin species have been identified, and they are referred to as cyclins A through I.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	GAK

<b>Immunogen:</b>	Synthesized peptide derived from Internal of human GAK
<b>Full name:</b>	cyclin G associated kinase
<b>Synonyms :</b>	DNAJ26; DNAJC26
<b>SwissProt:</b>	O14976
<b>WB Predicted band size:</b>	143 kDa
<b>WB Positive control:</b>	293 cells lysate
<b>WB Recommended dilution:</b>	500-3000

